
Controls – Ranger 3000


Button	Primary Function (Short Press <sup>1</sup> )	Secondary Function (Long Press <sup>2</sup> )	Menu Function (Short Press)
<div style="border: 1px solid black; padding: 2px; text-align: center;">On/Zero Off</div> <p style="text-align: center;">Yes</p>	<ul style="list-style-type: none"> <li>• Turns scale on</li> <li>• Sets display to zero</li> </ul>	<ul style="list-style-type: none"> <li>• Turns the scale off</li> </ul>	<ul style="list-style-type: none"> <li>• Selects/accepts displayed setting</li> </ul>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Print Units</div> <p style="text-align: center;">No</p>	<ul style="list-style-type: none"> <li>• Sends the current value to the selected COM ports if AUTOPRINT is set to Off.</li> </ul>	<ul style="list-style-type: none"> <li>• Toggles through active weighing units</li> </ul>	<ul style="list-style-type: none"> <li>• Advances to the next menu or menu item.</li> <li>• Rejects the current setting on the display and advances to the next available setting</li> </ul>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Function Mode</div> <p style="text-align: center;">Back</p>	<ul style="list-style-type: none"> <li>• Initiates an application mode</li> </ul>	<ul style="list-style-type: none"> <li>• Allows changing the application mode</li> </ul>	<ul style="list-style-type: none"> <li>• Returns to previous menu item</li> </ul>
<div style="border: 1px solid black; padding: 2px; text-align: center;">M+ Menu</div> <p style="text-align: center;">Exit</p>	<ul style="list-style-type: none"> <li>• Accumulates the weight and displays the accumulated information at 0 load.</li> </ul>	<ul style="list-style-type: none"> <li>• Enters User menu</li> </ul>	<ul style="list-style-type: none"> <li>• Quickly exit User Menu.</li> <li>• Aborts calibration in progress</li> </ul>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Tare</div>	<ul style="list-style-type: none"> <li>• Enter/clear a tare value. Clears the accumulation when the accumulation information is displayed.</li> </ul>		

**Notes:** <sup>1</sup> Short Press: Press less than 1 second.


<sup>2</sup> Long Press: Press and hold for more than 2 seconds.

## Controls – Ranger Count 3000

Button	Primary Function (Short Press)	Secondary Function (Long Press)	Menu Function (Short Press)	Library Function (Short Press)
 Yes	<ul style="list-style-type: none"> <li>• Turn scale on</li> <li>• Zero the display</li> </ul>	<ul style="list-style-type: none"> <li>• Turn scale off.</li> </ul>	<ul style="list-style-type: none"> <li>• Accept the current menu or setting.</li> </ul>	<ul style="list-style-type: none"> <li>• Accept the current setting.</li> </ul>
 No	<ul style="list-style-type: none"> <li>• Send the displayed value to the COM port</li> </ul>	<ul style="list-style-type: none"> <li>• Change the weighing unit.</li> </ul>	<ul style="list-style-type: none"> <li>• Advance to the next menu or setting.</li> <li>• Increment the displayed value.</li> </ul>	<ul style="list-style-type: none"> <li>• Advance to the next library or setting.</li> <li>• Increment the displayed value.</li> </ul>
 Back	<ul style="list-style-type: none"> <li>• Initiate the function of the current application mode.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch between Check Weigh, Check Count and Off.</li> </ul>	<ul style="list-style-type: none"> <li>• Previous menu</li> <li>• Decrement the displayed value.</li> </ul>	<ul style="list-style-type: none"> <li>• Previous library</li> <li>• Decrement the displayed value.</li> </ul>
 Exit	<ul style="list-style-type: none"> <li>• Accumulate the weight or pieces.</li> </ul>	<ul style="list-style-type: none"> <li>• Enter user menu.</li> <li>• View the Audit Trail event counters</li> </ul>	<ul style="list-style-type: none"> <li>• Exit the menu.</li> <li>• Abort the calibration</li> </ul>	<ul style="list-style-type: none"> <li>• Exit the library.</li> </ul>
	<ul style="list-style-type: none"> <li>• ID number input</li> </ul>			
	<ul style="list-style-type: none"> <li>• Display/Store APW</li> </ul>			
	<ul style="list-style-type: none"> <li>• Enter numeric values on the display.</li> </ul>			
 1 ΔΔ 2	<ul style="list-style-type: none"> <li>• Enter dot (.) on the display.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch between scale 1 and scale 2</li> </ul>		
 C +/-	<ul style="list-style-type: none"> <li>• Clear last character from the display.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch between positive and negative value</li> </ul>		
	<ul style="list-style-type: none"> <li>• Perform a tare operation.</li> </ul>			
	<ul style="list-style-type: none"> <li>• Display/Start APW</li> </ul>			

Controls – Valor 7000

Button	Primary Function (Short Press <sup>1</sup> )	Secondary Function (Long Press <sup>2</sup> )	Menu Function (Short Press)
 Yes	<ul style="list-style-type: none"> <li>• Turns scale on</li> <li>• Sets display to zero</li> </ul>	<ul style="list-style-type: none"> <li>• Turns the scale off</li> </ul>	<ul style="list-style-type: none"> <li>• Selects/accepts displayed setting</li> </ul>
 No	<ul style="list-style-type: none"> <li>• Sends the current value to the selected COM ports if AUTOPRINT is set to Off.</li> </ul>	<ul style="list-style-type: none"> <li>• Toggles through active weighing units</li> </ul>	<ul style="list-style-type: none"> <li>• Toggles through available settings</li> </ul>
 Back	<ul style="list-style-type: none"> <li>• Selects function setting</li> </ul>	<ul style="list-style-type: none"> <li>• Selects active Mode</li> </ul>	<ul style="list-style-type: none"> <li>• Returns to previous settings</li> </ul>
 Exit	<ul style="list-style-type: none"> <li>• Accumulates the weight and displays the accumulated information at 0 load.</li> </ul>	<ul style="list-style-type: none"> <li>• Enters User Menu</li> </ul>	<ul style="list-style-type: none"> <li>• Quickly exit User Menu</li> </ul>
	<ul style="list-style-type: none"> <li>• Enter/clear a tare value. Clears the accumulation when the accumulation information is displayed.</li> </ul>		
	<ul style="list-style-type: none"> <li>• The IR Sensor can be programmed to act as "touchless" button.</li> </ul>		

**Notes:** <sup>1</sup> Short Press: Press less than 1 second.  
<sup>2</sup> Long Press: Press and hold for more than 2 seconds.

Ohaus Corporation, 7 Campus Drive, Suite 310, Parsippany, New Jersey, 07054, USA www.ohaus.com

**Declaration of conformity** We, Ohaus Corporation, declare under our sole responsibility that the balance models listed below marked with "CE" – are in conformity with the directives and standards mentioned.

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
For non-automatic weighing instruments used in an Article 1, 2.(a) application, additional metrological marking according to Annex IV of Council directive 2009/23/EC must be attached to the instrument.

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Issued on: 2012-12-04

  
 Ted Xia  
 President  
 Ohaus Corporation  
 Parsippany, NJ USA

  
 Robert Hansen  
 Compliance Manager  
 Ohaus Corporation  
 Parsippany, NJ USA